

Product Name: Coverage Spray TB Plus

Product Code: 1629-08, 1629-B4

<b>INCI</b>	<b>CAS #</b>	<b>Purpose</b>
Aqua	7732-18-5	Water
Isopropyl Alcohol 99%, Anhydrous	67-63-0	Active Anti-Microbial
Butoxyethanol	111-76-2	Surfactant
Alkyl C12-18 dimethylbenzyl ammonium chloride	68391-01-5	Active Anti-Microbial
Alkyl C12-18 dimethyl ethylbenzyl ammonium chloride	85409-23-0	"
Sodium nitrite	7632-00-0	Corrosion inhibitor
Ethanol	64-17-5	Solvent
Sulfate of ethoxylated alcohol	Withheld	Cleaning agent
Nonionic surfactant	Withheld	Cleaning agent